




HSF1 (Ab-142) Antibody

Product Code	CSB-PA268961
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00613
Immunogen	Synthesized non-phosphopeptide derived from Human HSF1 around the phosphorylation site of threonine 142 (L-L-T(p)-D-V).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total HSF1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Heat shock factor protein 1; Heat shock transcription factor 1; HSF 1; HSTF 1; HSTF1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSF1

Image



Western blot analysis of extracts from K562 cells, using HSF1 (Ab-142) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.