




AKT1/AKT3 (Ab-437/434) Antibody

Product Code	CSB-PA932880
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31749
Immunogen	Synthesized non-phosphopeptide derived from Human AKT1/3 around the phosphorylation site of tyrosine 437/434 (T-R-Y(p)-F-D).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total AKT1/3 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AKT; AKT1 kinase; C-AKT; EC 2.7.11.1; kinase Akt1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKT1/AKT3

Image



Western blot analysis of extracts from A549 cells, using AKT1/3 (Ab-437/434) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.