





ICK (Ab-159) Antibody

Product Code	CSB-PA153658
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPZ9
Immunogen	Synthesized non-phosphopeptide derived from Human ICK around the phosphorylation site of tyrosine 159 (T-D-Y(p)-V-S).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total ICK protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.22; hICK; ICK; Intestinal cell kinase; KIAA0936
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ICK

Image



Western blot analysis of extracts from 3T3 cells, treated with PBS (10uM, 60mins), using ICK (Ab-159) antibody.



Western blot analysis of extracts from 293 cells (Lane 2) and cos-7 cells (Lane 3), using ICK (Ab-159) antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.