




INSR (Ab-1375) Antibody

Product Code	CSB-PA945193
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06213
Immunogen	Synthesized non-phosphopeptide derived from Human Stathmin around the phosphorylation site of threonine 1375 (I-L-T(p)-L-P).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total INSR protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CD220 antigen; EC 2.7.10.1; insulin receptor; IR; kinase InsR
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	INSR

Image



Western blot analysis of extracts from HepG2 cells, using INSR (Ab-1375) antibody.