




MAK (Ab-159) Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA056397 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P20794 |
| Immunogen | Synthesized peptide derived from Internal of Human MAK. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Specificity | The antibody detects endogenous levels of total MAK protein. |
| Tested Applications | ELISA,WB;WB:1:500-1:3000 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | EC 2.7.11.22; Male germ cell-associated kinase; Serine/threonine-protein kinase MAK; |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MAK |

Image



Western blot analysis of extracts from K562 cells, using MAK (Ab-159) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.