




MAP2K3 (Ab-222) Antibody

Product Code	CSB-PA084775
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46734
Immunogen	Synthesized non-phosphopeptide derived from Human MAP2K3 around the phosphorylation site of threonine 222 (A-K-T(p)-M-D).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total MAP2K3 protein
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	dual specificity mitogen-activated protein kinase kinase 3; EC 2.7.12.2; kinase MKK3; MAP kinase kinase 3; MAP2K3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAP2K3

Image



Western blot analysis of extracts from Jurkat cells, treated with serum (20%, 15mins), using MAP2K3 (Ab-222) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.