





RPS6KA5 (Ab-212) Antibody

Product Code	CSB-PA938251
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75582
Immunogen	Synthesized peptide derived from Internal of Human MSK1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MSK1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.1; Nuclear mitogen- and stress-activated protein kinase 1; 90 kDa ribosomal protein S6 kinase 5; RSK-like protein kinase; RSKL
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RPS6KA5

Image



Western blot analysis of extracts from mouse brain cells, using MSK1 (Ab-212) antibody.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using MSK1 (Ab-212) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.