




RPS6KA4 (Ab-568) Antibody

Product Code	CSB-PA903004
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75676
Immunogen	Synthesized non-phosphopeptide derived from Human MSK2 around the phosphorylation site of threonine 568 (M-Q-T(p)-P-C).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total MSK2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.1; kinase MSK2; KS6A4; Ribosomal protein kinase B; RPS6KA4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RPS6KA4

Image



Western blot analysis of extracts from RAW264.7 cells, treated with UV (5mins), using MSK2 (Ab-568) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.