




PAK2 (Ab-192) Antibody

Product Code	CSB-PA101227
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13177
Immunogen	Synthesized non-phosphopeptide derived from Human PAK2 around the phosphorylation site of serine 192 (T-K-S(p)-I-Y).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total PAK2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.1; gamma-PAK; kinase PAK2; P21-activated kinase 2; p21-activated kinase 2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PAK2

Image



Western blot analysis of extracts from HepG2 cells, treated with serum (20%, 15mins), using PAK2 (Ab-192) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.