





PAK3 (Ab-154) Antibody

Product Code	CSB-PA554964
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75914
Immunogen	Synthesized non-phosphopeptide derived from Human PAK3 around the phosphorylation site of serine 154 (Y-M-S(p)-F-T).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PAK3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	beta-PAK; CDC42/RAC effector kinase PAK-B; EC 2.7.1.; EC 2.7.11.1; kinase PAK3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PAK3

Image



Western blot analysis of extracts from rat heart, using PAK3 (Ab-154) antibody.



Western blot analysis of extracts from HeLa cells (Lane 2), using PAK3 (Ab-154) antibody. The lane on the left is treated with synthesized peptide.