




PRKCD (Ab-64) Antibody

Product Code	CSB-PA634286
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05655
Immunogen	Synthesized non-phosphopeptide derived from Human PKCD around the phosphorylation site of tyrosine 64 (H-I-Y(p)-E-G).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total PKCD protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.13; kinase PKC-delta; KPCD; nPKC-delta; PKC-delta
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKCD

Image



Western blot analysis of extracts from 3T3 cells, treated with EGF (200ng/ml, 30mins), using PKCD (Ab-64) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.