




PRKG2 (Ab-126) Antibody

Product Code	CSB-PA843898
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13237
Immunogen	Synthesized non-phosphopeptide derived from Human PKG2 around the phosphorylation site of serine 126 (G-V-S(p)-A-E).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PKG2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CGK 2; cGKII; cGMP-dependent protein kinase 2; EC 2.7.11.12; KGP2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKG2

Image



Western blot analysis of extracts from 293 cells, treated with anisomycin (25µg/ml, 30mins), using PKG2 (Ab-126) antibody.