




AR (Ab-363) Antibody

Product Code	CSB-PA082237
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10275
Immunogen	Synthesized non-phosphopeptide derived from Human Androgen Receptor around the phosphorylation site of tyrosine 363 (D-Y-Y(p)-N-F).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total Androgen Receptor protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ANDR; Androgen receptor; DHTR; Dihydrotestosterone receptor; NR3C4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AR

Image



Western blot analysis of extracts from K562 cells, treated with EGF (200ng/ml, 5mins), using Androgen Receptor (Ab-363) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.