




AR (Ab-363) Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA082237 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P10275 |
| Immunogen | Synthesized non-phosphopeptide derived from Human Androgen Receptor around the phosphorylation site of tyrosine 363 (D-Y-Y(p)-N-F). |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse |
| Specificity | The antibody detects endogenous levels of total Androgen Receptor protein. |
| Tested Applications | ELISA,WB;WB:1:500-1:3000 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | ANDR; Androgen receptor; DHTR; Dihydrotestosterone receptor; NR3C4 |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | AR |

Image



Western blot analysis of extracts from K562 cells, treated with EGF (200ng/ml, 5mins), using Androgen Receptor (Ab-363) antibody.