




PXN (Ab-272) Antibody

Product Code	CSB-PA287904
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49023
Immunogen	Synthesized non-phosphopeptide derived from Human Paxillin around the phosphorylation site of serine 272 (M-A-S(p)-L-S).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total Paxillin protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Paired box protein Pax-1; PAX-1; PAXI; paxillin; paxillin alpha
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PXN

Image



Western blot analysis of extracts from HeLa cells, using Paxillin (Ab-272) antibody.