




FOS (Ab-232) Antibody

Product Code	CSB-PA578446
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01100
Immunogen	Synthesized non-phosphopeptide derived from Human FOS around the phosphorylation site of threonine 232 (V-A-T(p)-P-E).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total FOS protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	activator protein 1; AP-1; C-FOS; cellular oncogene c-fos; Cellular oncogene fos
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOS

Image



Western blot analysis of extracts from RAW264.7 cells, using FOS (Ab-232) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.