




FGFR1/FGFR2 (Ab-463) Antibody

Product Code	CSB-PA160014
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11362
Immunogen	Synthesized peptide derived from Internal of Human FGFR1/2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total FGFR1/2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Basic fibroblast growth factor receptor 1; bFGF-R-1; BFGFR; CD331; CEK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FGFR1/FGFR2

Image



Western blot analysis of extracts from JurKat cells, using FGFR1/2 (Ab-463) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.