




BCL6 (Ab-333) Antibody

Product Code	CSB-PA935424
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41182
Immunogen	Synthesized non-phosphopeptide derived from Human Bcl-6 around the phosphorylation site of serine 333 (P-Q-S(p)-P-Q).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total Bcl-6 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	B-cell CLL/lymphoma 6; B-cell lymphoma 6 protein; B-cell lymphoma 6 protein homolog; Bcl-5; Bcl5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BCL6

Image



Western blot analysis of extracts from 293 cells, using Bcl-6 (Ab-333) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.