




BRAF (Ab-753) Antibody

Product Code	CSB-PA921529
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15056
Immunogen	Synthesized non-phosphopeptide derived from Human B-Raf around the phosphorylation site of threonine 735 (P-K-T(p)-P-I).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total B-Raf protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	B-RAF; B-RAF proto-oncogene serine/threonine-protein kinase; B-Raf proto-oncogene serine/threonine-protein kinase; BRAF; BRAF1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BRAF

Image



Western blot analysis of extracts from rat brain cells, using B-Raf (Ab-753) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.