





CLDN7 (Ab-210) Antibody

Product Code	CSB-PA046041
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95471
Immunogen	Synthesized non-phosphopeptide derived from Human Claudin 7 around the phosphorylation site of tyrosine 210 (S-K-E-Y(p)-V).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total Claudin 7 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CEPTRL2; claudin 7; claudin 9; claudin-1; claudin-7
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CLDN7


Image



Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue using Claudin 7 (Ab-210) antibody.



Western blot analysis of extracts from rat liver cells, using Claudin 7 (Ab-210) antibody.



Western blot analysis of extracts from 3T3 cells (Lane 2), using Claudin 7 (Ab-210) antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.