





MAPKAPK2 (Ab-272) Antibody

Product Code	CSB-PA876694
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49137
Immunogen	Synthesized non-phosphopeptide derived from Human MAPKAPK2 around the phosphorylation site of serine 272 (A-I-S(p)-P-G).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total MAPKAPK2 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.1; kinase MAPKAPK2; MAP kinase-activated protein kinase 2; MAPK-activated protein kinase 2; MAPK2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAPKAPK2

Image



Western blot analysis of extracts from COS cells, using MAPKAPK2 (Ab-272) antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using MAPKAPK2 (Ab-272) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.