




CNOT2 (Ab-101) Antibody

Product Code	CSB-PA951939
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZN8
Immunogen	Synthesized non-phosphopeptide derived from Human CNOT2 around the phosphorylation site of serine 101 (S-L-S(p)-Q-G).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total CNOT2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CCR4-associated factor 2; CCR4-NOT transcription complex subunit 2; CCR4-NOT transcription complex; subunit 2; CDC36
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CNOT2

Image



Western blot analysis of extracts from HepG2 cells, treated with starved (24hours), using CNOT2 (Ab-101) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.