




TAL1 (Ab-122) Antibody

Product Code	CSB-PA276394
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17542
Immunogen	Synthesized non-phosphopeptide derived from Human TAL-1 around the phosphorylation site of serine 122 (Q-L-S(p)-P-P).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total TAL-1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	T-cell acute lymphocytic leukemia protein 1; SCL; TCL5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TAL1

Image



Western blot analysis of extracts from Jurkat cells, treated with PMA (125ng/ml, 30mins), using TAL-1 (Ab-122) antibody.