




HCLS1 (Ab-397) Antibody

Product Code	CSB-PA790699
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14317
Immunogen	Synthesized peptide derived from C-terminal of Human HCLS1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total HCLS1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	cortactin-like; CTTNL; HCLS1; Hematopoietic cell-specific LYN substrate 1; Hematopoietic lineage cell-specific protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HCLS1

Image



Western blot analysis of extracts from JK cells (Lane 2), using HS1 (Ab-397) antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.