




ATF4 (Ab-219) Antibody

Product Code	CSB-PA123639
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18848
Immunogen	Synthesized peptide derived from C-terminal of Human ATF4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total ATF4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Activating transcription factor 4; activating transcription factor 4 (tax-responsive enhancer element B67); ATF4; cAMP response element-binding protein 2; cAMP-dependent transcription factor ATF-4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATF4

Image



Western blot analysis of extracts from A549 cells and Jurkat cells, using ATF-4 (Ab-219) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.