




MERTK/TYRO3 (Ab-753) Antibody

Product Code	CSB-PA114581
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12866
Immunogen	Synthesized peptide derived from Internal of Human MERTK.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total MERTK protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	c-mer; c-mer proto-oncogene tyrosine kinase; MER; MER receptor tyrosine kinase; MERTK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MERTK/TYRO3

Image



Western blot analysis of extracts from 3T3 cells (Lane 2), using MER (Ab-753) antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.