




ACTB Antibody

Product Code	CSB-PA284599
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60709
Immunogen	Synthesized peptide derived from Human Actin.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total Actin protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ACT; Actin; cytoplasmic 1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACTB

Image



Western blot analysis of extracts from HeLa cells, using Actin-pan antibody.



Immunohistochemical analysis of paraffin-embedded human skeletal muscle tissue using Actin-pan antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.