





APAF1 Antibody

Product Code	CSB-PA109273
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14727-6
Immunogen	Synthesized peptide derived from Human APAF-1-ALT.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total APAF-1-ALT protein.
Tested Applications	ELISA, WB, IHC, IF; WB:1:500-1:3000, IHC:1:50-1:200, IF:1:100-1:200
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	APAF1-interacting protein; APIP; APIP2; CGI-29 protein; MMRP19
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	APAF1

Image



Western blot analysis of extracts from COLO205, using APAF-1-ALT antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using APAF-1-ALT antibody.