





BAX Antibody

Product Code	CSB-PA167247
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07812
Immunogen	Synthesized peptide derived from Human Bax.
Raised In	Rabbit
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous levels of total Bax protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:200,IF:1:100-1:200
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Apoptosis regulator BAX; membrane isoform alpha; BAXA; Bcl2-associated X protein;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAX

Image



Western blot analysis of extracts from 293 cells, using Bax antibody.



Immunohistochemical analysis of paraffin-embedded human brain tissue using Bax antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.