





COL4A1 Antibody

Product Code	CSB-PA975386
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02462
Immunogen	Synthesized peptide derived from Human Collagen IV.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total Collagen IV protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	COL4A1 NC1 domain; COLLAGEN OF BASEMENT MEMBRANE; ALPHA-1 CHAIN ARRESTEN; alpha 1 type IV collagen preproprotein; collagen IV alpha-1 polypeptide
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COL4A1

Image



Western blot analysis of extracts from HeLa cells, using Collagen IV antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Collagen IV antibody.