





CCNA1 Antibody

Product Code	CSB-PA574179
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78396
Immunogen	Synthesized peptide derived from Human Cyclin A.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total Cyclin A protein.
Tested Applications	ELISA, WB, IHC, IF; WB: 1:500-1:3000, IHC: 1:50-1:100, IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cyclin A1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CCNA1


Image




Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using Cyclin A antibody.



Immunofluorescence analysis of COS7 cells, using Cyclin A antibody.



Western blot analysis of extracts from COS7 cells, using Cyclin A antibody.



Western blot analysis of extracts from JK cells (Lane 2) and cos-7 cells (Lane 3), using Cyclin A antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.