





TFDP1 Antibody

Product Code	CSB-PA149329
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14186
Immunogen	Synthesized peptide derived from Human DP-1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total DP-1 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	E2F DIMERIZATION PARTNER 1; TRANSCRIPTION FACTOR DP1; TRANSCRIPTION FACTOR; SEQUENCE-SPECIFIC; DRTF1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TFDP1


Image



Western blot analysis of extracts from HeLa cells, using DP-1 antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using DP-1 antibody.



Immunofluorescence analysis of HeLa cells, treated with Forskolin (40nM, 30mins), using DP-1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.