





E2F4 Antibody

Product Code	CSB-PA290646
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16254
Immunogen	Synthesized peptide derived from Human E2F4.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous levels of total E2F4 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	E2F-4; p107/p130-binding protein;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	E2F4


Image



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using E2F4 antibody.



Western blot analysis of extracts from COS7 cells, using E2F4 antibody.



Western blot analysis of extracts from JK cells (Lane 2), using E2F4 antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.