





E2F6 Antibody

Product Code	CSB-PA945555
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75461
Immunogen	Synthesized peptide derived from Human E2F6.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total E2F6 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	E2F-6; E2F transcription factor 6 isoform 1; E2F transcription factor 6 isoform 2; E2F transcription factor 6 isoform 3; E2F transcription factor 6 isoform 4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	E2F6


Image



Immunohistochemical analysis of paraffin-embedded human brain tissue using E2F6 antibody.



Western blot analysis of extracts from HepG2 cells, using E2F6 antibody.



Western blot analysis of extracts from Colo cells (Lane 2), using E2F6 antibody. The lane on the left is treated with synthesized peptide.