





FGFR1OP Antibody

Product Code	CSB-PA579407
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95684
Immunogen	Synthesized peptide derived from Human FGFR1 Oncogene Partner.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total FGFR1 Oncogene Partner protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	FOP; FIBROBLAST GROWTH FACTOR RECEPTOR 1 ONCOGENE PARTNER;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FGFR1OP

Image



Western blot analysis of extracts from HepG2 cells, using FGFR1 Oncogene Partner antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using FGFR1 Oncogene Partner antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.