





FHIT Antibody

Product Code	CSB-PA105645
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49789
Immunogen	Synthesized peptide derived from Human FHIT.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total FHIT protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AP3Aase; FRA3B; bis(5-adenosyl)-triphosphatase; inucleosidetriphosphatase; AP3A HYDROLASE FRAGILE SITE 3p14.2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FHIT

Image



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using FHIT antibody.



Immunofluorescence analysis of A549 cells, using FHIT antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.