





GRM7 Antibody

Product Code	CSB-PA962170
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14831
Immunogen	Synthesized peptide derived from Human GluR7.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total GluR7 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EAA5; GLR7; GluR7a; dJ1090M5.1 (glutamate receptor; ionotropic
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GRM7

Image



Western blot analysis of extracts from HuvEc cells, using GluR7 antibody.



Immunohistochemical analysis of paraffin-embedded human brain tissue using GluR7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.