




Phospho-INSR (Thr1375) Antibody

Product Code	CSB-PA994563
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06213
Immunogen	Peptide sequence around phosphorylation site of threonine 1375(I-L-T(p)-L-P) derived from Human IR .
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of INSR only when phosphorylated at threonine 1375.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	IR; CD220 antigen; kinase InsR;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	INSR

Image



Western blot analysis of extracts from HepG2 cells using INSR (Phospho-Thr1375) Antibody. The lane on the right is treated with the antigen-specific peptide.

Product Modify	Phospho-Thr1375
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.