





GUCY1B3 Antibody

Product Code	CSB-PA214476
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02153
Immunogen	Synthesized peptide derived from Human Guanylate Cyclase β .
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous levels of total Guanylate Cyclase β protein.
Tested Applications	ELISA, WB, IHC, IF; WB: 1:500-1:3000, IHC: 1:50-1:200, IF: 1:100-1:200
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	4.6.1.2; GC-S-beta-1; GC-SB3; GUC1B3; GUCB3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GUCY1B1

Image



Western blot analysis of extracts from HuvEc cells, using Guanylate Cyclase β antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using Guanylate Cyclase β antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.