





HDAC1 Antibody

Product Code	CSB-PA157026
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13547
Immunogen	Synthesized peptide derived from Human HDAC1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total HDAC1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Histone deacetylase 1; HD1; RPD3L1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC1

Image



Western blot analysis of extracts from NIH/3T3 cells, using HDAC1 antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using HDAC1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.