





HDAC7 Antibody

Product Code	CSB-PA094126
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUI4
Immunogen	Synthesized peptide derived from Human HDAC7.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total HDAC7 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	histone deacetylase 7A;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC7

Image



Western blot analysis of extracts from HuvEc cells, using HDAC7 antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using HDAC7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.