





DNAJB1 Antibody

Product Code	CSB-PA288820
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25685
Immunogen	Synthesized peptide derived from Human HSP40.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total HSP40 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	HSPF1; DNAJ E. COLI HOMOLOG OF SUBFAMILY B MEMBER 1; HEAT-SHOCK 40-KD PROTEIN 1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DNAJB1


Image



Western blot analysis of extracts from COLO205 cells, using HSP40 antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using HSP40 antibody.



Immunofluorescence analysis of NIH/3T3 cells, using HSP40 antibody.