





IVL Antibody

Product Code	CSB-PA796319
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07476
Immunogen	Synthesized peptide derived from Human Involucrin.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Involucrin protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Involucrin; IVL;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IVL

Image



Western blot analysis of extracts from HuvEc cells, using Involucrin antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Involucrin antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.