





MMP7 Antibody

Product Code	CSB-PA256113
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09237
Immunogen	Synthesized peptide derived from Human MMP-7.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total MMP-7 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	PUMP1; PUTATIVE METALLOPROTEINASE I; matrix metalloproteinase 7 (matrilysin; uterine);
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MMP7


Image



Immunohistochemical analysis of paraffin-embedded human colon carcinoma tissue using MMP-7 antibody.



Western blot analysis of extracts from COS7 cells, using MMP-7 antibody.



Western blot analysis of extracts from COS7 cells (Lane 2) and HuVc cells (Lane 3), using MMP-7 antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.