





MMP9 Antibody

Product Code	CSB-PA142806
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14780
Immunogen	Synthesized peptide derived from Human MMP-9.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MMP-9 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	92 kDa gelatinase; 92 kDa type IV collagenase precursor; CLG4B; GELB; Gelatinase B
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MMP9

Image



Western blot analysis of extracts from HeLa cells, using MMP-9 antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using MMP-9 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.