





NFE2L2 Antibody

Product Code	CSB-PA828758
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16236
Immunogen	Synthesized peptide derived from Human NRF2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total NRF2 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	HEBP1; NF-E2 related factor 2; NF2L2; NFE2-related factor 2; NFE2L2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NFE2L2


Image



Western blot analysis of extracts from HuvEc cells, using NRF2 antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using NRF2 antibody.



Immunofluorescence analysis of HuvEc cells, using NRF2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.