




PTH Antibody

Product Code	CSB-PA562344
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01270
Immunogen	Synthesized peptide derived from Human Parathyroid Hormone antibody.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Parathyroid Hormone protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	PTH/PTHr receptor; PTH/PTHrP type I receptor; PTHR1; Parathyroid hormone/parathyroid hormone-related peptide receptor precursor;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTH

Image



Western blot analysis of extracts from A549 cells, using Parathyroid Hormone antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.