





PTCH1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA914197 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q13635 |
| Immunogen | Synthesized peptide derived from Human Patched. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse |
| Specificity | The antibody detects endogenous levels of total Patched protein. |
| Tested Applications | ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | DROSOPHILA HOMOLOG OF 1; patched homolog (Drosophila); |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | PTCH1 |


Image



Immunohistochemical analysis of paraffin-embedded human brain tissue using Patched antibody.



Western blot analysis of extracts from mouse muscle cells, using Patched antibody.



Western blot analysis of extracts from 3T3 cells (Lane 2) and HepG2 cells (Lane 3), using Patched antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.