





PTCH1 Antibody

Product Code	CSB-PA914197
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13635
Immunogen	Synthesized peptide derived from Human Patched.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total Patched protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DROSOPHILA HOMOLOG OF 1; patched homolog (Drosophila);
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTCH1


Image



Immunohistochemical analysis of paraffin-embedded human brain tissue using Patched antibody.



Western blot analysis of extracts from mouse muscle cells, using Patched antibody.



Western blot analysis of extracts from 3T3 cells (Lane 2) and HepG2 cells (Lane 3), using Patched antibody. The lane on the left is treated with synthesized peptide.