




Phospho-PHB (Tyr259) Antibody

Product Code	CSB-PA242284
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35232
Immunogen	Peptide sequence around phosphorylation site of tyrosine 259 (I-T-Y(p)-L-P) derived from Human PHB.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of PHB only when phosphorylated at tyrosine 259.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using
Clonality	Polyclonal
Alias	Prohibitin; PHB1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHB

Image



Western blot analysis of extracts from MCF, PC12, 293 and 3T3 cells using PHB (Phospho-Tyr259) Antibody.

Product Modify Phospho-Tyr259